



The Aesthetic Lithium Disilicate System
Pressing • Layering • Staining



Shade Systems

LD PORCELAIN

SHADE	W1	W2	W3
OPAQUE DENTIN		-	OD-A1; OD-N=1:1
BODY		-	
OPAL ENAMEL	OPAL 56	OPAL 56: OPAL 57 =2:1	OPAL 56: OPAL 57 =1:1

A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
OD-A1	OD-A2	OD-A3	OD-A3.5	OD-A4	OD-B2; OD-N=1:1	OD-B2	OD-B2; OD-B4=1:1	OD-B4	OD-C2; OD-N=1:1	OD-C2	OD-C2; OD-C4=1:1	OD-C4	OD-D3; OD-N=2:1	OD-D3	OD-D3; OD-B4=1
A1B	A2B	A3B	A3.5B	A4B	B1B	B2B	B3B	B4B	C1B	C2B	C3B	C4B	D2B	D3B	D4B
OPAL 57	OPAL 58	OPAL 59	OPAL 60	OPAL 57	OPAL 58	OPAL 59	OPAL 60	OPAL 58	OPAL 59	OPAL 60	OPAL 57	OPAL 58	OPAL 59	OPAL 60	OPAL 57

OPAQUE DENTIN	OD-N														
OPAL ENAMEL	OPAL 56	OPAL 57	OPAL 58	OPAL 59	OPAL 60										
ENAMEL EFFECT	T	T-Glass	OPAL T	OPAL SL	BT	PT	GT	YT	OT	AM-Y					

GUM COLOR	GUM-1	GUM-2	GUM-3	GUM-4	GUM-5
CORRECTION	ADD-ON B	ADD-ON T			

* Available from September 2015

LD PRESS

SHADE	W1	W2	W3
VINTAGE LD PRESS		-	
	MT-W1	MT-W2	MT-W3

A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
T-1		T-2		T-3	T-1	T-2	T-3		T-1		T-3		T-1		T-3
MT-A1	MT-A2	MT-A3	MT-A3.5	MT-A4	MT-B1	MT-B2	MT-B3	MT-B4	MT-C1	MT-C2	MT-C3	MT-C4	MT-D2	MT-D3	MT-D4
MO-1		MO-2			MO-1		MO-2		MO-1		MO-3		MO-1		

ART LF STAIN

SHADE	IMAGE	EXAMPLE OF USE				
		SHADE ADJUSTMENT	CERVICAL	STAIN	INCISAL	GUM
P		●				●
Y		●	●			
Bl		●			●	
O		●	●	●		●
V		●			●	●
Gr		●		●		
B		●		●	●	
w					●	
DR-Br			●	●		
Br			●	●		
O-Br			●	●		
B-Br			●	●		

SHADE	IMAGE	EXAMPLE OF USE				
		SHADE ADJUSTMENT	CERVICAL	STAIN	INCISAL	GUM
K			●	●		
G		●			●	
Bl-G					●	
CY			●			●
RP			●			●
WR			●			
Vn					●	
AS		●	●			A1 A2 A3 A3.5 A4 rootA
BS		●	●			B1 B2 B3 B4 rootB
CS		●	●			C1 C2 C3 C4 rootC
DS		●	●			D2 D3 D4